




Howard Instruments, Inc.

68000-68425 DISPOSABLE TREPHINE BLADES SHORT - STERILE



- Available in 0.25mm increments, 6mm thru 9.5mm
- Available in 0.5mm increments, 2.5mm thru 5.5mm and 9.5mm thru 11mm
- Length 9mm, 0.350 inches

REGULAR SHORT TREPHINES

68000	6.00mm	68200	8.00mm
68025	6.25mm	68225	8.25mm
68050	6.50mm	68250	8.50mm
68075	6.75mm	68275	8.75mm
68100	7.00mm	68300	9.00mm
68125	7.25mm	68325	9.25mm
68150	7.50mm	68350	9.50mm
68175	7.75mm		

68600L-68975L DISPOSABLE TREPHINE BLADES LONG - STERILE



- Available in 0.25mm increments, 6mm thru 9.5mm
- Available in 0.5mm increments, 2.5mm thru 5.5mm and 10mm thru 16mm
- Length 16mm, 0.625 inches

REGULAR LONG TREPHINES

68600-L	6.00mm	68800-L	8.00mm
68625-L	6.25mm	68825-L	8.25mm
68650-L	6.50mm	68850-L	8.50mm
68675-L	6.75mm	68875-L	8.75mm
68700-L	7.00mm	68900-L	9.00mm
68725-L	7.25mm	68925-L	9.25mm
68750-L	7.50mm	68950-L	9.50mm
68775-L	7.75mm		



SPECIAL LONG TREPHINES

68525-L	2.50mm	68975-L	10.00mm	68998-L	13.50mm
68530-L	3.00mm	68980-L	10.50mm	69000-L	14.00mm
68535-L	3.50mm	68985-L	11.00mm	69002-L	14.50mm
68540-L	4.00mm	68987-L	11.50mm	69005-L	15.00mm
68545-L	4.50mm	68992-L	12.00mm	69007-L	15.50mm
68550-L	5.00mm	68994-L	12.50mm	69010-L	16.00mm
68555-L	5.50mm	68996-L	13.00mm		

TREPHINES/HANDLES

70000H HOWARD UNIVERSAL TREPHINE HANDLE



- For use with the series of short trephine blades
- By turning the back of the handle you can raise or lower the mm size that the handle accommodates
- With the center of the shaft made of solid metal, it allows a **POSITIVE LOCK** on the trephine. This eliminates movement of the trephine at the end of the handle.

Trephines from 6mm to 9.5mm are placed here.

Trephines from 10mm to 11mm are placed here.

70000BP BASE PLUG TO BE USED WITH CORNEAL PUNCH 22mm diameter X 9mm depth



70000BPS BASE PLUG SHIELD TO COVER CORNEAL TISSUE (FITS OVER BASE PLUG) 27mm diameter X 22mm depth



70000P HOWARD PUNCH

- To be used with Universal handles 70000H and 70016H
- Has self-aligning mechanism for easy insertion of the Howard Universal Trephine Blade



Website: howardinstruments.com
 email: howard@howardinstruments.com
 Phone: (205) 553-4453 Fax: (205) 56 9267

