



INFUSION CANNULAS

DISPOSABLE INFUSION CANNULA

- 20 Gauge
- 45 degree beveled tip
- Luer lock adapter
- Uses 260mm of Silastic® tubing
- Length 270mm, 10.6 inches

IS-1025

- 2.5mm

IS-1040

- 4.0mm

IS-1060

- 6.0mm



METAL INFUSION CANNULAS

HIGH VISCOSITY DISPOSABLE METAL INFUSION CANNULA

- 20 Gauge
- 45 degree beveled tip
- Can be used as standard infusion cannula
- Designed for easy injection of high viscosity substance such as silicone oil
- Tubing length 381mm, 15 inches

MIS-1040

- 4.0mm tip

MIS-1060

- 6.0mm tip



SILICONE OIL / PERFLOURDECALINE

SILICONE OIL

NOTFDAAPPROVED

SO-1000

- Sterile
- 1000 centistoke
- 10ml

SO-2000

- Sterile
- 2000 centistoke
- 10ml

SO-5000

- Sterile
- 5000 centistoke
- 10ml



PERFLOURDECALINE

NOTFDAAPPROVED

PFD-1000

- Sterile
- 5ml

PFD-2000

- Sterile
- 7ml


Howard Instruments, Inc.
 4749 Appletree
 Tuscaloosa, AL 35405
 Website: howardinstruments.com
 email: howard@howardinstruments.com
 Phone: (205) 553-4453 Fax: (205) 556 9267



Howard Instruments, Inc. *SPECIAL VITRO-RETINAL ACCESSORIES*



H-2800 Viscous Fluid Injector



H-2860 Perflourdecaline Injection
Needle

SILICONE OIL / PERFLOURDECALINE

SILICONE OIL NOT FDA APPROVED

SO-1000

- Sterile
- 1000 centistoke
- 10ml

SO-2000

- 2000 centistoke

SO-5000

- 5000 centistoke



PERFLOURDECALINE NOT FDA APPROVED

PFD-1000

- Sterile
- 5ml

PFD-2000


- 7ml




REUSABLE SCLERAL PLUGS & TABLET



R20-10000
20 gauge



R19-10000
19 gauge




R18-10000
18 gauge

Supplied with 4 plugs
Color-coded
Steam Autoclavable silicone rubber tablet
Raised support for easy grip with forceps (not included)

DISPOSABLE SCLERAL PLUGS & TABLET

RDP-1000

Supplied with 4 plugs,
2, 19 gauge
2, 20 gauge
Polypropylene tablet



RCA-1000 Scleral Plug Forceps

- Lightweight
- Developed for use with virtually any traditional plug set available
- Manufactured from the highest quality American stainless steel