

IrisMate[®] Iris Retractors

Flexible Translimbal Iris Retractors



“A Superior Product at a Fraction of the Cost”

Iris Retractor 2000

CATALOG NUMBER IRIS-2000

- Each Disposable Case contains 4 disposable Iris Retractors
- Specially designed orange cinch for easier viewing and superior grip.
- Bent end allows for proper positioning under microscope.
- Easy Insertion and Removal
- All MEDICAL GRADE materials
- Flexible nylon hook retracts iris with minimal damage to the tissue
- Design for working through 0.3mm - 0.5mm self sealing incision
- Safely used in Aphakic, Pseudophakic and Phakic Eye



CALL or EMAIL TODAY to Place Your Order!