

Infusion Cannula

Disposable Infusion Cannula (Sterile)

- IS-1025** 2.5mm Beveled tip, packaged in box of ten
IS-1040 4.0mm Beveled tip, packaged in box of ten
IS-1060 6.0mm Beveled tip, packaged in box of ten
- 20 Gauge
 - Luer lock adapter
 - Uses 260mm of Silastic® tubing
 - Length 270mm, 10.6 inches

Disposable Infusion Cannula (Non Sterile)

- IS-1025-NS** 2.5mm Beveled tip, packaged in packages of 100
IS-1040-NS 4.0mm Beveled tip, packaged in packages of 100
IS-1060-NS 6.0mm Beveled tip, packaged in packages of 100
- 20 Gauge
 - Luer lock adapter
 - Uses 260mm of Silastic® tubing
 - Length 270mm, 10.6 inches



High Viscosity Disposable Infusion Cannula(Sterile)

MIS-1040

- 4.0mm tip, packaged in box of ten
- 20 Gauge 45 degree beveled tip
- Can be used as standard infusion cannula.
- Designed for easy injection of high viscosity substance.
- Tubing length 10 inches

High Viscosity Disposable Infusion Cannula (Sterile)

MIS-1060-6

- 6.0mm tip, packaged in box of ten
- 20 Gauge 45 degree beveled tip
- Can be used as standard infusion cannula.
- Designed for easy injection of high viscosity substance.
- Tubing length 6 inches

High Viscosity Disposable Infusion Cannula(Non Sterile)

MIS-1060-NS

- 6.0mm tip, packaged in packages of 100
- 20 Gauge 45 degree beveled tip
- Can be used as standard infusion cannula.
- Designed for easy injection of high viscosity substance.

High Viscosity Disposable Infusion Cannula(Sterile)

MIS-1040-6

- 4.0mm tip, packaged in box of ten
- 20 Gauge 45 degree beveled tip
- Can be used as standard infusion cannula.
- Designed for easy injection of high viscosity substance.
- Tubing length 6 inches

High Viscosity Disposable Infusion Cannula (Sterile)

MIS-1060

- 6.0mm tip, packaged in box of ten
- 20 Gauge 45 degree beveled tip
- Can be used as standard infusion cannula.
- Designed for easy injection of high viscosity substance.
- Tubing length 10 inches

High Viscosity Disposable Infusion Cannula(Non Sterile)

MIS-1040-NS

- 4.0mm tip, packaged in packages of 100
- 20 Gauge 45 degree beveled tip
- Can be used as standard infusion cannula.
- Designed for easy injection of high viscosity substance.
- Tubing length 100mm