


  
**HOWARD INSTRUMENTS II, INC.**

---

# *IrisMate<sup>®</sup> Iris Retractors*

Flexible Translimbal Iris Retractors




## **Iris Retractor 1000**

**CATALOG NUMBER**    **IRIS-1000**

Each Sterilization Case contains 6 Reusable Iris Retractors

- Affordable
- Easy Insertion and Removal
- All MEDICAL GRADE materials
- Flexible nylon hook retracts iris with minimal damage to the tissue
- Specially designed orange silicone cinch provides superior grip
- Works through 0.3mm - 0.5mm self sealing incision
- Safely used in Aphakic, Pseudophakic and Phakic Eyes



## **Iris Retractor 2000**

**CATALOG NUMBER**    **IRIS-2000**

Each Disposable Case contains 4 disposable Iris Retractors

- Cinch is color coded orange for easier viewing under microscope.
- Bent end allows for proper positioning under microscope.
- Affordable
- Easy Insertion and Removal
- All MEDICAL GRADE materials
- Flexible nylon hook retracts iris with minimal damage to the tissue
- Specially designed orange silicone cinch provides superior grip
- Works through 0.3mm - 0.5mm self sealing incision
- Safely used in Aphakic, Pseudophakic and Phakic Eye

