

HOWARD INSTRUMENTS II, INC.

DISPOSABLE TREPHINE BLADES / LONG - STERILE (68600L-68970L)

Long Trephines available in 0.25mm increments, 6mm thru 9.75mm
 Special Long Trephines available in 0.5mm increments,
 2.5mm thru 5.5mm and 10mm thru 16mm
 Length 16mm, 0.625 inches



LONG TREPHINES

68600-L	6.00mm	68800-L	8.00mm
68625-L	6.25mm	68825-L	8.25mm
68650-L	6.50mm	68850-L	8.50mm
68675-L	6.75mm	68875-L	8.75mm
68700-L	7.00mm	68900-L	9.00mm
68725-L	7.25mm	68925-L	9.25mm
68750-L	7.50mm	68950-L	9.50mm
68775-L	7.75mm	68970-L	9.75mm

SPECIAL LONG TREPHINES

68525-L	2.50mm	68975-L	10.00mm	68998-L	13.50mm
68530-L	3.00mm	68980-L	10.50mm	69000-L	14.00mm
68535-L	3.50mm	68985-L	11.00mm	69002-L	14.50mm
68540-L	4.00mm	68987-L	11.50mm	69005-L	15.00mm
68545-L	4.50mm	68992-L	12.00mm	69007-L	15.50mm
68550-L	5.00mm	68994-L	12.50mm	69010-L	16.00mm
68555-L	5.50mm	68996-L	13.00mm		

69030-L 18.00mm diameter x 16.00mm length , Designed for Eye Banks with a Special Knurled Handle

DISPOSABLE TREPHINE BLADES / SHORT - STERILE (68000-68350)

Available in 0.25mm increments and 6mm thru 9.75mm
 Length 9mm, 0.350 inches




SHORT TREPHINES

68000	6.00mm	68200	8.00mm
68025	6.25mm	68225	8.25mm
68050	6.50mm	68250	8.50mm
68075	6.75mm	68275	8.75mm
68100	7.00mm	68300	9.00mm
68125	7.25mm	68325	9.25mm
68150	7.50mm	68350	9.50mm
68175	7.75mm	68375	9.75mm

TREPHINE ACCESSORIES


70000H
HOWARD UNIVERSAL TREPHINE HANDLE

- For use with the 68000 series trephine blades
- By turning the back of the handle you can raise or lower the mm size that the handle accommodates
- With the center of the shaft made of solid metal, it allows a **POSITIVE LOCK** on the trephine. This eliminates movement of the trephine at the end of the handle.



70000BP
BASE PLUG
TO BE USED WITH CORNEAL PUNCH
22mm dia. X 9mm depth

70000BPS
BASE PLUG SHIELD
TO COVER CORNEAL TISSUE
(FITS OVER BASE PLUG)
27mm dia. X 22mm depth



70000P
HOWARD PUNCH
TO BE USED WITH CORNEAL HANDLES
70000H and 70016H

